

**List of thermoplastic disks D. 134 mm
Thermoforming Temperature and Time
for BB93 equipment**

Gruppo	Spess.	Tipo	Codice	Temperatura	Tempo
Cappette	0,5	Calcinable	190002	1.9 (200°C)	30"
Base Plate	1,5	Base Plate	190305	1.5 (175°C)	1,15"
	2		190306	1.5 (175°C)	1,40"
Porta impronte	3	Tray	190013	1.9 (200°C)	3'
	4		190014	1.9 (200°C)	4'
Provvi-sori	5	Tray	190015	1.9 (200°C)	5'
	0,5		190340	1.5 (175°C)	25"
Provvi-sori neg.	1	Temporary A/A	190342	1.5 (175°C)	55"
	1,5		190325	1.5 (175°C)	1,15"
Provvi-sori neg.	2	Temporary negative	190326	1.5 (175°C)	1.40"
	1		190329	1.3 (150°C)	45"
Byte Plane Morbidi	2	Byte Soft	190330	1.3 (150°C)	1,30"
	3		190331	1.3 (150°C)	2,20"
	4		190332	1.3 (150°C)	3'
	0,5	Byte Hard	190347	1.5 (175°C)	25"
Byte Plane Rigi	0,75		190346	1.5 (175°C)	40"
	1		190341	1.5 (175°C)	55"
	2		190335	1.5 (175°C)	1,40"
	3		190336	1.5 (175°C)	2,20"
	4		190337	1.5 (175°C)	3'
Sbiancamento	0,75	Bleaching Hard	190364	1.5 (175°C)	40"
	1		190365	1.5 (175°C)	55"
	1	Bleaching Soft	190368	1.3 (150°C)	45"
	0,75	Bleach.SHS(S25+H25+S25)	190051	1.9 (200°C)	45"
Byte Plane R-M Dopo 30 sec di riscaldamento riserrare l'oblo	2	Byte HS (H15 + S5)	190052	1.9 (200°C)	2'
	3	Byte HS (H20 + S10)	190053	1.9 (200°C)	3'
	4	Byte HS (H20 + S20)	190054	1.9 (200°C)	4'
	1,5	Muonthguard	190071	1.3 (150°C)	1,15"
Paradenti Sportivi	2		190072	1.3 (150°C)	1,30"
	3		190074	1.3 (150°C)	2,20"
	4		190075	1.3 (150°C)	3'