

- Sandblaster -

"SANDY"

SANDY is a sandblaster specially designed and built to meet the needs of the lab attached to the dental office.

• ERGONOMIC • FUNCTIONAL • COMPACT

SANDY presents a serie of innovations studied to optimize its use and to remedy to daily problems. Ergonomically adapted to the normal working position of a man. The operator has the control of all the commands which are easy to operate; the appliance is compact but it provides a big internal working area. The sandblaster must be connected to the aspirator AS800 to remove the dust formed during the processing. The joint of connection to the aspirator can be placed both to the right and to the left allowing excellent placement in constricted spaces.



It displays the operating pressure.

PRESSURE REGULATOR

It activates the operating pressure from 0,5 to 6 bar. Attached is the filter for moisture and impurities.

ON - OFF

It activates the functioning of the solenoid and the inside light.

HOLDERS

Standard: 25÷50μm 50÷125μm On demand: 125÷250μm

ASPIRATOR CONNECTION

Diameter 32 mm can be placed both to the right and to the left.

SELECTOR

To select the type of sand to use.

PEDAL POWER

Through the solenoid it drives the output of the sand from the chosen noozle.

ANTI-SCRATCH LID

The lid of the working area is made with a special anti-scratch glass.

NOOZLES

Specific pointed to work in small cavities.

